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
1. Opentrons Labworks Inc. 2. seqWell, Inc.

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

- ExpressPlex™ Library Prep Kit Single-Step Workflow Demo for High-Throughput Multiplexed Sample Prep. (2023). YouTube. <https://www.youtube.com/watch?v=gB4IsVH0AIA>
- Leonard, J. T., Feeley, R., Huang, S., Wang, L., & Smith, J. (n.d.). A transformational one-step library preparation method for multiplexed plasmid and amplicon sequencing. Library Prep Simplified through Scalable Workflows | seqWell. https://seqwell.com/wp-content/uploads/2022/06/20220519_ExpressPlex-Poster.pdf
- Leonard, J. T., Huang, S., Feeley, R., Couture, J., Palys, J., Smith, J., & Liu, Y. (n.d.). An Optimal Strategy for Internalizing Ultra-high Throughput Plasmid Sequencing. Library Prep Simplified through Scalable Workflows | seqWell. https://seqwell.com/wp-content/uploads/2023/05/SW2023_Optimal-Strategy-HTP-Sequencing-11x8.5-05192023.pdf
- Palys, J., Pearson, A., Whitmore, I., & Rees, H. (2023). High-Throughput NGS Library Prep for Plasmids and Amplicons with the ExpressPlex™ Library Preparation Kit on the SPT Labtech firefly®.
- SeqWell Inc. (2023). ExpressPlex™ Library Prep Kit. Library Prep Simplified through Scalable Workflows | seqWell. https://seqwell.com/wp-content/uploads/2023/04/SW2023_ExpressPlex-Product-Flyer-04202023.pdf